

**CHEM / BCMB 4190/6190/8189**


**Introductory NMR**

Lecture 13

## The DEPT Experiment

### Interpreting $^{13}\text{C}$ NMR spectra:

- Very useful to know which signals belong to quaternary carbons, CH, CH<sub>2</sub>, and CH<sub>3</sub>
- J-modulated spin-echo experiment and refocused INEPT are useful but don't provide all the necessary information
- The DEPT experiment (Distorsionless Enhancement by Polarization Transfer) is now one of the most important techniques available for interpreting 1D  $^{13}\text{C}$  NMR spectra



**Figure 8-22.**

The DEPT experiment; for pulse sequence see text.

A: 100.6 MHz  $^{13}\text{C}$  NMR spectrum of the neuraminic acid derivative **1** with  $^1\text{H}$  BB decoupling ( $\delta = 10$  to 110 region only).

B: CH sub-spectrum: DEPT(90).

C: CH<sub>2</sub> sub-spectrum: DEPT(45) – DEPT(135).

D: CH<sub>3</sub> sub-spectrum: DEPT(45) + DEPT(135) – 0.707 DEPT(90).

(Experimental conditions:

167 mg of the compound in 2.3 ml

D<sub>2</sub>O; 10 mm sample tube;

16 K data points; 32 FIDs for

$\theta_1 = 45^\circ$  and  $\theta_3 = 135^\circ$ , 64 FIDs for

$\theta_2 = 90^\circ$ ;  $\tau = 3.57$  ms; total time

approx. 12 min.)

## Pulse Sequence


$^1\text{H}$  channel:

$90^\circ_{x'} - \tau - 180^\circ_{x'} - \tau - \Theta_{y'} - \tau - \text{BB decoupling}$

$^{13}\text{C}$  channel:

$90^\circ_{x'} - \tau - 180^\circ - \tau - \text{FID } (t_2)$

- $\tau$  is  $1/2J_{\text{CH}}$  with  $J_{\text{CH}} = 140 \text{ Hz}$
- is chosen to be  $\Theta = 45^\circ$ ,  $\Theta = 90^\circ$ ,  $\Theta = 135^\circ$
- Vector diagrams are not adequate to explain the effect of this pulse sequence and to understand this pulse sequence.
- The following curves explain the effect of  $\Theta$  on the intensities of signals




**Figure 8-23.**

DEPT experiment. Curves calculated from Equation (8-6) for the intensities of CH, CH<sub>2</sub> and CH<sub>3</sub> signals as functions of the pulse angle  $\Theta_{y'}$ ; CH: - · - · - · - · -, CH<sub>2</sub>: ———, CH<sub>3</sub>: - - - - -.

- CH sub-spectrum:  $\Theta_2 = 90^\circ$  or DEPT (90)
- CH<sub>2</sub> sub-spectrum: DEPT(45) - DEPT(135)
- CH<sub>3</sub> sub-spectrum: DEPT(45) + DEPT(135) - 0.707DEPT(90)
- To compare absolute intensities the DEPT(90) must be collected with twice as many scans as for the other two DEPT experiments

Practically, it is sufficient to carry out two experiments, DEPT(90) and DEPT(135)



**Figure 8-24.**

DEPT(135) spectrum of the neuraminic acid derivative **1**, recorded using the pulse sequence given in the text, with  $\theta_y = 135^\circ$ . The signals of the five CH groups, identified with the help of the DEPT(90) spectrum ( $\theta_y = 90^\circ$ ), are marked by arrows. The other three positive signals arise from CH<sub>3</sub> groups, and the two negative signals from CH<sub>2</sub> groups.

(Experimental conditions:

20 mg of the compound in 0.5 ml D<sub>2</sub>O; 5 mm sample tube; 32 K data points; 300 FIDs;  $\tau = 3.57$  ms; total time approx. 20 min.)

**Table 8-2.**

Partial assignment of the <sup>13</sup>C NMR signals of **1** from the results of the DEPT experiment.

$\delta$	CH <sub>3</sub>	CH <sub>2</sub>	CH	C	Assignment
23.2	×				CH <sub>3</sub> (Ac)
40.31		×			C-3
52.12	×				
52.83			×		C-5
54.65	×				
64.50		×			C-9
67.51			×		
69.18			×		
70.98			×		
71.67			×		
100.32				×	C-2
171.50 <sup>a)</sup>				×	
175.93 <sup>a)</sup>				×	

<sup>a)</sup> Values from the complete spectrum (Fig. 8-12 A).