


RELAXATION

REVIEW: TRANSVERSE (SPIN-SPIN, T_2) RELAXATION

- Our oscillating signal decays as a function of time as the phase coherence between the precessing magnetic dipoles (vectors) is lost
- Transverse (spin-spin) or T_2 relaxation refers to loss of phase coherence in the transverse (x - y) plane
- This loss of coherence is due to local magnetic field differences experienced by the nuclei







- T_2 relaxation is characterized both by a component due to magnetic field inhomogeneity (the uninteresting component) and a by component arising from local fluctuating magnetic fields produced by the nuclei themselves

$$\frac{1}{T_2^*} = \frac{1}{T_2} + \frac{1}{T_{2(\Delta B_0)}}$$

a


REVIEW: TRANSVERSE (T_2) RELAXATION

**b****c****d**

- The decay of signal due to T_2 relaxation is first order




$$\frac{dM_y}{dt} = -\frac{M_y}{T_2^*}$$

$$M_y = M_{y_0} e^{(-t/T_2^*)}$$



REVIEW: TRANSVERSE (T_2) RELAXATION

- Short T_2 times lead to broad lines (undesirable, lower signal-to-noise)
- A poorly shimmed magnet leads to short T_2 times and broad peaks




- The Fourier transform of a decaying exponential function ($\cos \omega t \cdot e^{(-t/T_2)}$) gives a **Lorentzian** line shape
- The width of a Lorentzian peak at 1/2 of the maximum height is $1/(\pi T_2^*)$

$$\Delta v_{1/2} = \frac{1}{\pi T_2^*} = \frac{1}{\pi T_2} + \gamma \Delta B_0$$

REVIEW: LONGITUDINAL (SPIN-LATTICE, T_1) RELAXATION

- Longitudinal (spin-lattice) or T_1 relaxation refers to return to thermal equilibrium (along the z axis) of the spin populations (return to equilibrium N_α and N_β values) following perturbation
- T_1 relaxation represents a loss of energy (heat) from the spins to the surroundings (i.e. the “lattice”: neighboring molecules, solvent, impurities, NMR tube walls, etc.)




- The return to equilibrium is first order
 - after a 90° pulse, the return to equilibrium is described as shown below:

$$\frac{dM_z}{dt} = -\frac{M_z - M_0}{T_1} \quad M_z = M_0(1 - e^{(-t/T_1)})$$

-after a 180° pulse.....

$$M_z = M_0(1 - 2e^{(-t/T_1)})$$

REVIEW: T_1 RELAXATION AND DATA ACQUISITION



- Need to wait for a time (" d_1 ") after acquisition of the FID in order for the spin populations to return to thermal equilibrium
 - after a 90° pulse, should wait for $\sim 5T_1$
 - thus, are limited in the number of "scans" that can be acquired in a given amount of time
- Important point: signal-to-noise (S/N) increases in proportion to the square root of the number of scans
$$S/N \propto \sqrt{N_S}$$
- For instance, if 1 scan gives a particular S/N, in order to double the S/N, 4 scans are required (unfortunately)

REVIEW: OPTIMIZING PULSE WIDTH

What is the optimal pulse width/length/angle to use?

- a 90° pulse angle gives maximum S/N for one pulse, but the delay (d_1) between successive pulses must be long for recovery of equilibrium magnetization
- following a 90° pulse, if $(d_1 + AQ) = 5 \times T_1$, then > 99% of equilibrium magnetization will be recovered before the next pulse
- thus, *for a given number of pulses* (without regards to time), 90° pulses will give maximum sensitivity (as long as $(d_1 + AQ) \geq 5 \times T_1$)
- other schemes, which permit faster pulsing (shorter d_1) combined with smaller pulse angles are possible
- for optimizing the S/N and total experimental time, the best compromise for the pulse width/angle is the *Ernst Angle*
- Ernst Angle is in degrees



$$\cos \alpha_{\text{Ernst}} = e^{-((d_1 + AQ)/T_1)}$$

REVIEW: OPTIMIZING PULSE WIDTH

$$\cos \alpha_{\text{Ernst}} = e^{-((d_1 + AQ)/T_1)}$$

What is the optimal repetition time (best S/N) for a given repetition time?

- the optimal pulse angle is dependent on the repetition time
- as the repetition time shortens, so does the optimal pulse angle




What is the optimal repetition time (best S/N) for a given experiment time using 90 degree pulses?

- the optimal pulse angle is $1.3 \times T_1$ for 90 degree pulses for a pre-determined total experiment time


MEASURING T_1 RELAXATION

- Longitudinal (spin-lattice) or T_1 relaxation refers to return to thermal equilibrium (along the z axis) of the spin populations (return to equilibrium N_α and N_β values) following perturbation
- The equilibrium values of N_α and N_β are perturbed following a pulse
 - following a 90° pulse, $N_\alpha = N_\beta$, and N_α and N_β return to equilibrium values with time



MEASURING T_1 RELAXATION


- If M_z could be measured with time following a 90° pulse, T_1 could be determined easily



- The return to equilibrium is first order
 - after a 90° pulse, the return to equilibrium is described as shown below:

$$\frac{dM_z}{dt} = -\frac{M_z - M_0}{T_1}$$

$$M_z = M_0(1 - e^{(-t/T_1)})$$



MEASURING T_1 RELAXATION

- If M_z could be measured with time following a 180° pulse, T_1 could be determined easily


$$M_z = M_0$$

$$M_z = -M_0$$

$$M_z \rightarrow M_0$$

$$M_z = M_0$$


$$\begin{aligned} N_\alpha &= N_{\alpha \text{ eq}} \\ N_\beta &= N_{\beta \text{ eq}} \\ N_{\alpha \text{ eq}} &> N_{\beta \text{ eq}} \end{aligned}$$



$$\begin{aligned} N_\alpha &= N_{\beta \text{ eq}} \\ N_\beta &= N_{\alpha \text{ eq}} \\ N_\beta &= N_{\alpha \text{ eq}} \end{aligned}$$


$$\begin{aligned} N_\alpha &\rightarrow N_{\alpha \text{ eq}} \\ N_\beta &\rightarrow N_{\beta \text{ eq}} \end{aligned}$$

$$\begin{aligned} N_\alpha &= N_{\alpha \text{ eq}} \\ N_\beta &= N_{\beta \text{ eq}} \\ N_{\alpha \text{ eq}} &> N_{\beta \text{ eq}} \end{aligned}$$




$$M_z = M_0(1 - 2e^{(-t/T_1)})$$

↑
note !





MEASURING T_1 RELAXATION: INVERSION RECOVERY

- Problem: **can't measure M_z directly**. Can only measure bulk magnetization in transverse plane
- Solution: ***inversion recovery experiment***
- The 180° pulse places the bulk magnetization vector along the $-z$ axis (NO components in the transverse plane)
- During the τ delay period, as the α and β spin populations return towards their equilibrium values (i.e. T_1 relaxation), the absolute magnitude of the bulk magnetization vector along $-z$ decreases initially (short τ values), becomes zero, and then increases along $+z$ (long τ values) until it equals M_0 (return to thermal equilibrium)
- The 90° “read” pulse places the magnetization in the transverse ($x-y$) plane, and thus the magnitude of the bulk magnetization vector in the transverse plane is proportional to the magnitude of the bulk vector when it was along z
- ***The magnitude of the Fourier transformed signal produced, therefore, is proportional to the magnitude of the bulk vector when it was along z***
- So, the intensity of the signal as a function of τ , where I_0 is the intensity when $\tau = \infty$, allows determination of T_1




$$M_z = M_0(1 - 2e^{(-t/T_1)}) \quad I = I_0(1 - 2e^{(-t/T_1)}) \quad I = I_0(1 - 2e^{(-\tau/T_1)})$$

INVERSION RECOVERY





here, the magnitude of bulk magnetization vector in transverse plane is proportional to the magnitude of the bulk magnetization vector along $\pm z$ after the τ period/delay

so, **signal intensity** after Fourier transform now is proportional to the magnitude of the bulk magnetization vector along $\pm z$ after the τ period/delay




MEASURING T_1 RELAXATION: INVERSION RECOVERY

- For well resolved signals, can determine all T_1 values in a single inversion recovery experiment (i.e. spectra acquired at many different τ values)
- For each signal, data (I versus “ τ ” or “ τ' ”) are fit to equation shown to get T_1 values





$$I = I_0(1 - 2e^{(-t/T_1)})$$




MEASURING T_1 RELAXATION: INVERSION RECOVERY

- Measuring T_1 for ^{13}C is no different
- Broadband ^1H -decoupling used usually to improve signal-to-noise



$$I = I_0(1 - 2e^{(-t/T_1)})$$



MEASURING T_1 RELAXATION: INVERSION RECOVERY

- Measure I versus time (t) or “ τ ”, fit to the equation
- Can make an *estimate* without fitting by determining τ when the signal intensity is 0 (i.e. when $I = 0$)

$$0 = I_0(1 - 2e^{(-t/T_1)}) = I_0 - I_0 2e^{(-t/T_1)}$$

$$I_0 = I_0 2e^{(-t/T_1)} \quad 1 = 2e^{(-t/T_1)} \quad 1/2 = e^{(-t/T_1)}$$

$$\ln(1/2) = -t/T_1 \quad \ln(2) = t/T_1 \quad T_1 = \frac{t}{\ln(2)}$$

usually written $T_1 = \frac{\tau_{\text{zero}}}{\ln(2)}$



- Typically, it is not advisable to linearize the equation
- However, the equation can be linearized and T_1 determined from the slope of the appropriate plot

$$I = I_0(1 - 2e^{(-t/T_1)}) = I_0 - I_0 2e^{(-t/T_1)}$$

$$I_0 - I = I_0 2e^{(-t/T_1)}$$


$$\ln(I_0 - I) = \ln 2I_0 - t/T_1 \quad (\text{slope} = -1/T_1)$$

$$I = I_0(1 - 2e^{(-t/T_1)})$$




MEASURING T_2 RELAXATION

- A 90° ($\pi/2$) pulse equilibrates N_α and N_β and promotes “**phase coherence**” among the nuclear magnetic dipoles, creating ***transverse*** magnetization.




- Transverse (spin-spin) or T_2 relaxation refers to the time-dependent loss of phase coherence in the transverse (x - y) plane



MEASURING T_2 RELAXATION


- T_2 relaxation is characterized both by a component due to magnetic field inhomogeneity (the *uninteresting* component) and a by component arising from local fluctuating magnetic fields produced by the nuclei themselves (the *interesting* component)
- The combination of these components is called T_2^*
- T_2^* is EASY to measure.....

$$\frac{1}{T_2^*} = \frac{1}{T_2} + \frac{1}{T_{2(\Delta B_0)}}$$




$$\frac{dM_y}{dt} = -\frac{M_y}{T_2^*}$$

$$M_y = M_{y_0} e^{(-t/T_2^*)}$$




MEASURING T_2^* RELAXATION

- T_2^* is EASY to estimate also
- Estimate M_y/M_{y0} at a particular time
 - estimate M_y/M_{y0} is approximately 0.5 at 3 seconds, so





- T_2^* can also be estimated from width of the Lorentzian NMR signal at $\frac{1}{2}$ of its maximum height

$$\Delta v_{1/2} = \frac{1}{\pi T_2^*}$$




MEASURING T_2 RELAXATION: SPIN-ECHO EXPERIMENT

- The natural component of T_2^* contributed by the local fluctuating magnetic fields of the nuclei themselves - the interesting component that we'll call T_2 - is more difficult to measure than T_2^*
- The contribution from magnetic field inhomogeneities has to be eliminated
- The most common way to do this is with a **spin-echo experiment**
- The spin echo experiment for measuring T_2 is based on the fact that the effects of magnetic field inhomogeneity can be **refocused**, and thus eliminated, leaving only the effects of the natural T_2 component to lead to loss of phase coherence




MEASURING T_2 RELAXATION: SPIN-ECHO EXPERIMENT

- In order to understand how the spin echo method works, we have to first consider how the 180° pulses affect magnetic dipole vectors in the transverse plane
- Following a 90° pulse, the bulk magnetization vector is in the x - y plane
- In the rotating frame, we'll assume that the Larmor frequency of the nucleus is slightly faster than the reference frequency, so, after a short time τ , the vector has moved away from the y axis and towards the x axis





MEASURING T_2 RELAXATION: SPIN-ECHO EXPERIMENT



- So, the net result in all cases is simply ***reflection through the z-x plane for a $180^\circ x$ pulse and reflection through the z-y plane for a $180^\circ y$ pulse***
- So, this is a simple way to decide the effect of a 180° pulse without considering the effect on each of the component vectors



MEASURING T_2 RELAXATION: SPIN-ECHO EXPERIMENT

- In the rotating frame, in the first τ period, the individual dipoles lose phase coherence due to T_2^* processes
 - The $180^\circ x$ pulse reflects the vectors through the z - x plane
 - In the next τ period, the dipole vectors refoc



MEASURING T_2 RELAXATION: SPIN-ECHO EXPERIMENT

- So, a spin-echo is an increased signal intensity resulting from refocusing of nuclear spin dipoles




- Between each pair of 180° pulses the dipoles refocus, producing a signal intensity increase or “echo”
- The 180° pulses refocus the effects of magnetic field inhomogeneity, but not true T_2 spin-spin relaxation processes


- T_2^*
- T_2
- spin “echos”
- Thus, each successive echo has a reduced intensity due to true T_2 relaxation processes

MEASURING T_2 RELAXATION: SPIN-ECHO EXPERIMENT

- The intensities of the echos depend only on T_2 , and not T_2^* , as a simple exponential



$$I = I_0 e^{(-t/T_2)}$$




- As with T_2^* , T_2 can be estimated by estimating I/I_0 at any given time

$$I = I_0 e^{(-t/T_2)} \quad \frac{I}{I_0} = e^{(-t/T_2)} \quad \ln\left(\frac{I}{I_0}\right) = \frac{-t}{T_2} \quad T_2 = -t / \ln\left(\frac{I}{I_0}\right)$$

- The equation can also be linearized, but this is best used for rough estimates, and should not replace fitting to the exponential equation


$$I = I_0 e^{(-t/T_2)} \quad \ln I = \ln(I_0) - t/T_2$$



MORE ON RELAXATION.....

T_1 – INFLUENCE OF DIRECTLY ATTACHED HYDROGENS


- The main contribution to spin-lattice (T_1) relaxation is dipole-dipole (through space) interactions
- For ^{13}C , the primary mechanism for relaxation is via the dipolar interaction with hydrogens
- Hydrogens attached directly to the ^{13}C nucleus contribute most significantly due to the close proximity between the hydrogens and ^{13}C nucleus
- T_1 values for ^{13}C nuclei decrease with increasing numbers of attached hydrogens



- The T_1 (${}^{-13}\text{CH}$) $\approx 2 T_1$ (${}^{-13}\text{CH}_2$) not perfect, but pretty close
- For $-\text{CH}_3$, fast rotation about the C-C bond leads to smaller τ_c and relatively larger T_1
- These depend on the magnetic moments (so, ^2H less efficient at promoting relaxation than ^1H)


T_1 – INFLUENCE OF NEIGHBORING HYDROGENS

- Through space, dipolar interactions promote relaxation
- These include interactions between nuclei not directly bonded to one another



3A

- The magnetic moment of ^2H (D) is less than the magnetic moment of ^1H
- Thus, ^2H (D) is much less efficient at promoting T_1 relaxation than ^1H




3B

T_1 – INFLUENCE OF MOLECULAR SIZE

- The rotational correlation time for a molecule, τ_c , is proportional to the inverse of the T_1 relaxation time constant

$$\tau_c \propto T_1^{-1}$$


polystyrene: molecular weight 100,000 – 400,000



- Relaxation times for large molecules are short (large τ_c)
- Relaxation times for small molecules are long (small τ_c)

T_1 – SEGMENTAL MOBILITIES




- Different parts of molecules have different mobilities, and, therefore, different effective correlation times
- This leads to variations in relaxation times due to these motional variations



- In straight-chain polymers, for example, the central regions exhibit relatively less motional flexibility, larger effective correlation times, and shorter T_1 relaxation times
- The ends of such polymers exhibit relative more motional flexibility, smaller effective correlation times, and longer T_1 relaxation times

T_1 – MOTIONAL ANISOTROPY

- Molecular reorientation or rotation does not necessarily occur isotropically (i.e. the same or uniformly in all directions or about all axes)
- Anisotropic motion can lead to distinct differences in correlation times for nuclei within a molecule, and distinctly different T_1 relaxation times



- In some benzene derivatives, correlation times for nuclei on the axis of rotation may be much longer than those for groups not on the axis
- Thus, nuclei on the axis of rotation display shorter T_1 relaxation times