


2D NMR Spectroscopy


Two-Dimensional NMR

- One-dimensional (1D) NMR: one frequency axis (x), and one intensity axis (y)
- Two-dimensional (2D) NMR: usually two frequency axes (x and y), and one intensity axis (z)
 - often, both frequency axes are chemical shifts (chemical shift correlation experiments)
 - one frequency axis, however, might be, for instance, a coupling constant (in Hz)
 - the directly detected dimension is f_2 , and the indirectly detected dimension f_1
- In two dimensional experiments, the x and y dimensions are related or *correlated* based on:
 - through bond coupling (COSY, TOCSY, HSQC)
 - through space interactions (NOESY)
 - chemical exchange processes (EXSY)




Two-Dimensional NMR


- One-dimensional experiments can be described as consisting of a spin “preparation” (P) phase (often just a single pulse) and a “detection” (D) phase, when the magnetization is detected




- Two-dimensional experiments include, in addition, a magnetization “evolution” period (E) and a “mixing” (M) period
 - during the evolution period, the modulation of the signal for the second dimension occurs
 - during the mixing period, magnetization or “coherence” transfer between spins occurs




- Additional E and M periods can be added to create 3D and 4D experiments



Some Two-Dimensional NMR Experiments




Two-Dimensional NMR Data Collection and Manipulation



- In a two-dimensional experiment, t_1 is incremented
 - a series of 1D experiments are collected at increasing t_1 values
 - these 1D FIDs are Fourier transformed (the “ t_2 ” transform) to give a series of 1D NMR spectra

Two-Dimensional NMR Data Collection and Manipulation


- Peaks (signals) in the series of 1D spectra are modulated by t_1 and decay with t_1
 - thus, the signals in the 1D spectra form a FID in the t_1 (second) dimension
 - so, a Fourier transform (the “ t_1 ” transform) can be performed on these signals to give the 2D plot
 - different signals in f_2 may be modulated differently in t_1 , resulting in different f_1 frequencies for different f_2 signals



COSY: 2D Chemical Shift Correlation



- COSY correlates chemical shifts of coupled nuclei
 - magnetization is transferred between coupled nuclei as a result of the COSY experimental scheme
 - provides information on through-bond connectivity

- consider two coupled spins, A and X:
 - X magnetization will precess at ν_X during t_1
 - X magnetization that is *not* subsequently transferred to A will also precess at ν_X during t_2
 - however, X magnetization subsequently transferred to A will precess during t_2 at ν_A
 - thus, off diagonal signals (with different chemical shifts in f_2 and f_1) indicate that X and A are coupled
 - diagonal elements are, in general, not interesting





COSY: Example

- COSY spectrum of an azo-sugar
- crosspeaks appear at positions corresponding to the chemical shifts of ^1H nuclei on adjacent carbon atoms
 - crosspeaks indicate significant coupling between the ^1H nuclei
 - *usually*, therefore, crosspeaks are not observed for ^1H nuclei separated by more than 3 bonds



TOCSY: “Total” Correlation Spectroscopy

- TOCSY spectrum of an azo-sugar
- The TOCSY “spin-lock” sequence enables transfer of magnetization between nuclei even if the coupling between them is very weak
 - thus, for the example shown, crosspeaks between most of the nuclei are observed



Heteronuclear Correlation Spectroscopy

- HSQC spectrum of menthol
- No diagonal peaks
- The chemical shifts of directly bonded ^1H - ^{13}C pairs are correlated
 - one axis is ^{13}C chemical shift
 - one axis is ^1H chemical shift

